FerroOrange

货号: F374 规格: 1 tube / 3 tubes

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

1. Maintaining Iron Homeostasis Is the Key Role of Lysosomal Acidity for Cell Proliferation, <i>Molecular Cell</i> , 2020 , 77(3), 645-655.e7	293T cells (人肾上皮细胞系)
2. Versatile iron-vitamin K3 derivative-based nanoscale coordination polymer augments tumor ferroptotic therapy, <i>Nano Research</i> , 2020 , doi.org/10.1007/s12274-020-3241-7	CT26 cells (小鼠结肠癌细胞)
3. Induction of ferroptosis in response to graphene quantum dots through mitochondrial oxidative stress in microglia, <i>Part. Fibre. Toxicol.</i> , 2020 , 17(1): 30. doi: 10.1186/s12989-020-00363-1	BV2 cells (小鼠小胶质瘤细胞)
4. FTY720 induces ferroptosis and autophagy via PP2A/AMPK pathway in multiple myeloma cells, <i>Life Sciences</i> , 2020 , 260, 118077	Multiple Myeloma cells (多发性骨髓瘤细胞)
5. PM2.5 induces ferroptosis in human endothelial cells through iron overload and redox imbalance, <i>Environmental Pollution</i> , 2019 , doi.org/10.1016/j.envpol.2019.07.105	Endothelial cells (血管内皮细胞)
6. Activation of p62-Keap1-Nrf2 Pathway Protects 6-Hydroxydopamine-Induced Ferroptosis in Dopaminergic Cells, <i>Molecular Neurobiology</i> , 2020 , 57(11), 4628-4641	SH-SY5Y cells (人神经母细胞瘤细胞)
7. Inhibition of ferroptosis protects House Ear Institute-Organ of Corti 1 cells and cochlear hair cells from cisplatin-induced ototoxicity, <i>J. Cell Mol. Med.</i> , 2020 , 24:12065-12081	HEI-OC1 cells (耳蜗毛细胞)
8. Cytochrome promotes sensitivity to ferroptosis by regulating p53 - YAP1 axis in colon cancer cells, <i>J. Cell. Mol. Med.</i> , 2021 , 00, 1-12	HCT 116 cells (人结肠癌细胞)
9. Dihydroartemisinin-loaded Magnetic Nanoparticles for Enhanced Chemodynamic Therapy, <i>Front. Pharmacol.</i> , 2020 , doi.org/10.3389/fphar.2020.00226	MCF-7/ADR cells (人乳腺癌细胞)
10. Periodontitis-level butyrate-induced ferroptosis in periodontal ligament fibroblasts by activation of ferritinophagy, <i>Cell Death Discovery</i> , 2020 , 6,119, doi.org/10.1038/s41420-020-00356-1	Periodontal Ligament Fibroblasts cells (牙周膜成纤维细胞)

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